

Durability. Security. **Control.**

Sharps container
cabinets designed to
minimize risk and
maximize security



TOUCHLESS DROP METHOD

Compared to standard mechanical door levers instructing users to “lift to assure disposal,” our touchless vertical drop method utilizes a drop tube that mitigates the possibility of touch contamination and cross contamination risks tied to high touch surfaces

DURABLE

Steel material (16 or 18 gauge) is durable, easy to clean, long lasting, and secure compared to any standard plastic enclosure

EASY TO MOUNT

Rugged back plate cut-outs are easy to mount to any wall surface yet secure to discourage forced removal

ANTI-BACTERIAL

Material options include 304 stainless steel, known as the industry standard for mitigating bacteria formation on high touch devices





CABINET COLORS & MATERIALS

Choose from a number of decorative colors to compliment a variety of interior environments or choose attractive polished stainless steel for indoor or outdoor placement

MATERIAL OPTIONS



*Stainless Steel
Carbonless*



*Stainless Steel
Polished*



Powder Coated

COLOR OPTIONS



Dark



Pearl



*Pediatric
Lavender*



CABINET LOCK OPTIONS

Post Medical utilizes solid steel material (not a metal veneer) secured behind a variety of lock options that include electronic key/locks (coming soon), tubular, or flat based on the security needs of the customer.



TUBULAR CAM LOCK

Tubular Key locks are ideal for customers who prefer a lower cost option with higher security access control



FLAT KEY

Flat key locks offer the most economical option offered to customers that are less concerned about access control measures



CYBERLOCK

Industry leading electronic lock technology offers access control capabilities ideal for securing Post Medical's sharps cabinets, medical field cases, pharmaceutical safes and more



1.5 GAL SHARPS CONTAINER CABINET

Description	SKU
Powder Coated, View Window	88-SFP
Powder Coated, No View Window	88-SFP-B
Powder Coated, Narrow View Window	88-SFP-C
Stainless Steel, View Cutout	88-SFP-SS
Stainless Steel All-Weather, View Cutout	88-SFP-SS-AW



1 QT SHARPS CONTAINER CABINET

Description	SKU
Stainless Steel	110-SS
Stainless Steel All-Weather	110-SS-AW
Powder Coated	110-S



1.5 QT SHARPS CONTAINER CABINET

Description	SKU
Stainless Steel	150-SS
Stainless Steel All-Weather	150-SS-AW
Powder Coated	150-SC



1-1.5 GAL SHARPS CONTAINER CABINET

Description	SKU
Stainless Steel	700-SS
Stainless Steel All-Weather	700-SS-AW
Powder Coated	700-A



2.5 GAL SHARPS CONTAINER CABINET

Description	SKU
Stainless Steel	750-SS
Stainless Steel All-Weather	750-SS-AW
Powder Coated	750-L