

Durability. Security. **Control.**

Sharps tabletop
container covers
designed to
minimize accidental
tipping and spillage



TOUCHLESS DROP METHOD

Compared to standard mechanical door levers instructing users to “lift to assure disposal,” our touchless vertical drop method utilizes a drop tube that mitigates the possibility of touch contamination and cross contamination risks tied to high touch surfaces

DURABLE

Steel material (16 or 18 gauge) is durable, easy to clean, long lasting, and secure compared to any standard plastic enclosure

EASY TO MOUNT

Rugged back plate cut-outs are easy to mount to any wall surface yet secure to discourage forced removal

ANTI-BACTERIAL

Material options include 304 stainless steel, known as the industry standard for mitigating bacteria formation on high touch devices

TABLE TOP COVERS

PRODUCT CODE DESCRIPTION QTY

CC-1.0-SS	CONSEAL COVER, FOR WD-110 CONTAINER, POWDER COAT STAINLESS STEEL FINISH	EACH
CC-1.0-W	CONSEAL COVER FOR WD-110 CONTAINER, POWDER COAT WHITE FINISH	EACH
CC-1.5-SS	CONSEAL COVER FOR WD-150 CONTAINER, POWDER COAT STAINLESS STEEL FINISH	EACH
CC-1.5-W	CONSEAL COVER FOR WD-150 CONTAINER, POWDER COAT WHITE FINISH	EACH
CC-2.0-SS	CONSEAL COVER FOR WD-200 CONTAINER, POWSER COAT STAINLESS STEEL FINISH	EACH
CC-2.0-W	CONSEAL COVER FOR WD-200 CONTAINER, POWDER COAT WHITE FINISH	EACH





1.5 GAL SHARPS CONTAINER CABINET

Description	SKU
Powder Coated, View Window	88-SFP
Powder Coated, No View Window	88-SFP-B
Powder Coated, Narrow View Window	88-SFP-C
Stainless Steel, View Cutout	88-SFP-SS
Stainless Steel All-Weather, View Cutout	88-SFP-SS-AW



1 QT SHARPS CONTAINER CABINET

Description	SKU
Steel 16 Gauge Amber	110-S
Stainless Steel	110-SST



1.5 QT SHARPS CONTAINER CABINET

Description	SKU
Steel 16 Gauge Amber	150-SC
Stainless Steel	150-SST



1.5 GAL SHARPS CONTAINER CABINET

Description	SKU
Steel Glacier White	700-DF-IT-W
Steel Grey	700-DF-IT-SS
Stainless Steel	700-DF-IT-SST



2.5 GAL SHARPS CONTAINER CABINET

Description	SKU
Steel Glacier White	750-IT-DF-W
Steel Grey	750-IT-DF-SS
Stainless Steel	750-IT-DF-SST