

# Durability. Security. Control.

Sharps tabletop container covers designed to minimize accidental tipping and spillage



## TOUCHLESS DROP METHOD

Post Medical Covers are engineered with a screw in bracket that secures directly onto the container's top for a seamless virtual drop of the used sharp directly into the container.

## DURABLE

Powder-Coated Steel material (18 gauge) is durable, easy to clean, long lasting, and secure compared to any standard plastic enclosure and long-lasting sharps accessory.

## STABLE

OSHA guidelines require "containers for disposable sharps to be kept upright to keep sharps and any liquids from spilling out of the container." The U.S. Department of Health and Human Services (HSS) officials recommend "the use of mechanisms to restrain sharps containers is one way of preventing spillage during use.

## ANTI-BACTERIAL

Powder-Coated Glacier White or Grey color options are easy to clean mitigating bacteria formation on high touch devices.

## DECORATIVE

Enhance the aesthetics of any environment with our attractive table top cover.

# PRODUCT SIZE, DIMENSIONS AND SPECIFICATIONS

| SKU       | Size        | Color         | Quantity | Length (Inches) | Width (Inches) | Height (Inches) | Compatible Container |
|-----------|-------------|---------------|----------|-----------------|----------------|-----------------|----------------------|
| CC-1.0-SS | 1qt Cover   | Grey Steel    | 1        | 5.8             | 6.5            | 6.5             | WD-110               |
| CC-1.0-W  | 1qt Cover   | Glacier White | 1        | 5.8             | 6.5            | 6.5             | WD-110               |
| CC-1.5-SS | 1.5qt Cover | Grey Steel    | 1        | 5.8             | 6.5            | 8               | WD-150               |
| CC-1.5-W  | 1.5qt Cover | Glacier White | 1        | 5.8             | 6.5            | 8               | WD-150               |
| CC-2.0-SS | 2qt Cover   | Grey Steel    | 1        | 5.8             | 6.5            | 10.5            | WD-200               |
| CC-2.0-W  | 2qt Cover   | Glacier White | 1        | 5.8             | 6.5            | 10.5            | WD-200               |

Patented

Made in the USA

 **CONSEAL!**

CONSEAL! CONTAINER COVERS



**1 QT COVER**

| Description   | SKU       |
|---------------|-----------|
| Grey Steel    | CC-1.0-SS |
| Glacier White | CC-1.0-W  |

**1.5 QT COVER**

| Description   | SKU       |
|---------------|-----------|
| Grey Steel    | CC-1.5-SS |
| Glacier White | CC-1.5-W  |



**2 QT COVER**

| Description   | SKU       |
|---------------|-----------|
| Grey Steel    | CC-2.0-SS |
| Glacier White | CC-2.0-W  |